La evidencia la
tiene...

 **smith&nephew**
TRIGEN[®]
INTERTAN[®]
Clavo anterógrado intertrocanterico

Menor

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Referencias:

1. American Academy of Orthopaedic Surgeons. Hip fractures in seniors: a call for health system reform. Position Statement 1144. Rosemont, IL: 1999.
2. Baumgaertner M, Solberg B. Awareness of tip-apex distance reduces failure of fixation of trochanteric fractures of the hip. *J Bone Joint Surg Br*. 1997;79:969-971;
3. Bentler SE, Liu L, Obrizan M, Cook EA, Wright KB, Geweke JF, et al. The aftermath of hip fracture: discharge placement, functional status change, and mortality. *Am. J Epidemiol*. 2009 Nov 15;170(10):1290-9. doi: 10.1093/aje/kwp266.
4. Erez O, Dougherty PJ. Early complications associated with cephalomedullary nail for intertrochanteric hip fractures. *The journal of trauma and acute care surgery*. Feb 2012;72(2):E101-105.
5. Galli M, Ciriello V, Bocchino L, Gangemi NM, Peruzzi M, Marzetti E. Clinical and functional outcomes of internal fixation with intertrochanteric antegrade nail in older patients with proximal extracapsular femoral fractures. *Eur J Trauma Emerg Surg*. 2013/10/17 2013:1-6.
6. Hoffmann S, Paetzold R, Stephan D, Püschel K, Buehren V, Augat P. Biomechanical evaluation of interlocking lag screw design in intramedullary nailing of unstable pertrochanteric fractures. *J Orthop Trauma*. 2013;27(9):483-490.
7. Hsueh KK, Fang CK, Chen CM, et al. Risk factors in cutout of sliding hip screw in intertrochanteric fractures: an evaluation of 937 patients. *Int Orthop*. 2010;34:1273-1276;
8. Kim JW, Kim TY, Ha YC, Lee YK, Koo KH. Outcome of intertrochanteric fractures treated by intramedullary nail with two integrated lag screws: A study in Asian population. *Indian J Orthop*. Jul-Aug 2015;49(4):436-441.
9. Liu Y, Tao R, Liu F, Wang Y, Zhou Z, Cao Y, et al. Mid-term outcomes after intramedullary fixation of peritrochanteric femoral fractures using the new proximal femoral nail antirotation (PFNA). *Injury* 2010;41:810-7. [http:// dx.doi.org/10.1016/j.injury.2010.03.020](http://dx.doi.org/10.1016/j.injury.2010.03.020).
10. Ma KL et al. Proximal femoral nails antirotation, Gamma nails, and dynamic hip screws for fixation of intertrochanteric fractures of femur: A Meta-Analysis. 2014. *Orthopaedics & traumatology: Surgery and Research* 100, 859-866.
11. Magaziner J, Fredman L, Hawkes W, Hebel JR, Zimmerman S, Orwig DL, Wehren L. Changes in functional status attributable to hip fracture: a comparison of hip fracture patients to community-dwelling aged. *Am J Epidemiol*. 2003; 157:1023-31.
12. Matre K, Vinje T, Havelin LI, et al. TRIGEN INTERTAN intramedullary nail versus sliding hip screw: a prospective, randomized multicenter study on pain, function, and complications in 684 patients with an intertrochanteric or subtrochanteric fracture and one year of follow-up. *The Journal of bone and joint surgery. American volume*. Feb 6 2013;95(3):200-208.
13. Mir HR, Edwards P, Sanders R, Haidukewych G. Iatrogenic displacement of minimally or nondisplaced intertrochanteric fractures during intramedullary fixation: incidence, potential etiologies, and clinical impact. *Current Orthopaedic Practice*. 2013;24(1):58-63.
14. Mir HR, Edwards P, Sanders R, Haidukewych G. Results of cephalomedullary nail fixation for displaced intracapsular femoral neck fractures. *Journal of orthopaedic trauma*. 2011;25(12):714-720.
15. Mundi S et al. Similar mortality rates in hip fracture patients over the past 31 years: A systematic review of RCTs. *Acta Orthopaedica* 2014; 85(1): 54-59.
16. Nayak A, Smithson I, Cooper S, Cox J, Marberry, S, Santoni BG, et al. Comparison of femoral head rotation and varus collapse between a single and integrated dual screw intertrochanteric hip fracture fixation device using a chair rise biomechanical model. Abstract presented at: Orthopaedic Trauma Association Annual Meeting; October 15-18, 2014; Tampa, FL, USA.
17. Ruecker AH, Rupprecht M, Gruber M, Gebauer M, Barvencik F, Briem D, et al. The treatment of intertrochanteric fractures: results using an intramedullary nail with integrated cephalocervical screws and linear compression. *J Orthop Trauma* 2009;23:22-30.
18. Rueger J, Moore C. Shortening of the femoral neck following peritrochanteric fracture. *Bone Joint Sci*. 2011 May;2(5).
19. Rupprecht M, Grossterlinden L, Ruecker AH, et al. A comparative biomechanical analysis of fixation devices for unstable femoral neck fractures: the Intertan versus cannulated screws or a dynamic hip screw. *J Trauma*. 2011;71:625-634
20. Sahin EK, Imerci A, Kinik H, Karapinar L, Canbek U, Savran A. Comparison of proximal femoral nail antirotation (PFNA) with AO dynamic condylar screws (DCS) for the treatment for unstable peritrochanteric femoral fractures. *Eur J Orthop Surg Traumatol* 2014;24:347-52.
21. Serrano-Riera R, Blair JA, Downes K, Sanders R. Cephalo-medullary nail fixation of intertrochanteric fractures: are two proximal screws better than one? Abstract presented at: Orthopaedic Trauma Association Annual Meeting; October 15-18, 2014; Tampa, FL, USA.
22. Seyhan M, Turkmen I, Unay K, Ozkut AT. Do PFNA devices and Intertan nails both have the same effects in the treatment of trochanteric fractures? A prospective clinical study. *J Orthop Sci*. Nov 2015;20(6):1053-1061.
23. Tao R, Lu Y, Xu H, Zhou ZY, Wang YH, Liu F. Internal fixation of intertrochanteric hip fractures: a clinical comparison of two implant designs. *Sci World J* 2013;2013:1-8.
24. Wang Q, Yang X, He HZ, Dong LJ, Huang DG. Comparative study of InterTAN and Dynamic Hip Screw in treatment of femoral intertrochanteric injury and wound. *International journal of clinical and experimental medicine*. 2014;7(12):5578-5582.
25. Wu D, Ren G, Peng C, Zheng X, Mao F, Zhang Y. InterTan nail versus Gamma3 nail for intramedullary nailing of unstable trochanteric fractures. *Diagnostic pathology*. 2014;9:191
26. Wu Y, Watson JT, Kuldjanov D, Jackman J. Rotationally stable fixation for intertrochanteric hip fractures: the Intertan experience, surgical technique, and outcomes. *Techniques in Ortho*. 2014; 29;3:120-132
27. Yu J et al. Internal fixation treatments for intertrochanteric fracture: a systematic review and meta-analysis of randomized evidence. *Nature Scientific Reports*. 2015.
28. Zhang S, Zhang K, Jia Y, Yu B, Feng W. InterTan nail versus Proximal Femoral Nail Antirotation-Asia in the treatment of unstable trochanteric fractures. *Orthopedics*. Mar 2013;36(3):e288-294.

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©Marcas registradas de Smith & Nephew.
Smith & Nephew, Inc.
©2016 Smith & Nephew, Inc.